


Secaubrytriol

Chemical Properties

CAS No. :	925932-10-1
Formula:	C ₃₀ H ₅₀ O ₅
Molecular Weight:	490.71
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Secaubrytriol is a natural product for research related to life sciences. The catalog number is TN5936 and the CAS number is 925932-10-1.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0379 mL	10.1893 mL	20.3786 mL
5 mM	0.4076 mL	2.0379 mL	4.0757 mL
10 mM	0.2038 mL	1.0189 mL	2.0379 mL
50 mM	0.0408 mL	0.2038 mL	0.4076 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481