

## Baogongteng C

## Chemical Properties


CAS No. :

Formula: C<sub>7</sub>H<sub>13</sub>NO<sub>2</sub>

Molecular Weight: 143.186

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Baogongteng C is a useful organic compound for research related to life sciences and the catalog number is T126450.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.9837 mL	34.9186 mL	69.8373 mL
5 mM	1.3967 mL	6.9837 mL	13.9675 mL
10 mM	0.6984 mL	3.4919 mL	6.9837 mL
50 mM	0.1397 mL	0.6984 mL	1.3967 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481