


## Antiviral agent 24

## Chemical Properties

CAS No. :	23661-03-2
Formula:	C <sub>18</sub> H <sub>18</sub> F <sub>3</sub> N <sub>5</sub> O <sub>4</sub>
Molecular Weight:	425.36
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Nucleoside Derivatives - 6-Modified purine nucleosides; Drugs and Inhibitors
Targets(IC50)	Others,Virus Protease

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3509 mL	11.7547 mL	23.5095 mL
5 mM	0.4702 mL	2.3509 mL	4.7019 mL
10 mM	0.2351 mL	1.1755 mL	2.3509 mL
50 mM	0.047 mL	0.2351 mL	0.4702 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481