


2',3',5'-Tri-O-benzoyl-4-thiouridine

Chemical Properties

CAS No. : 15049-50-0
 Formula: C₃₀H₂₄N₂O₈S
 Molecular Weight: 572.59
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	2',3',5'-Tri-O-benzoyl-4-thiouridine is a Nucleoside Derivative - Thio-nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7465 mL	8.7323 mL	17.4645 mL
5 mM	0.3493 mL	1.7465 mL	3.4929 mL
10 mM	0.1746 mL	0.8732 mL	1.7465 mL
50 mM	0.0349 mL	0.1746 mL	0.3493 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481