


Luviquat mono CP

Chemical Properties

CAS No. : 85563-48-0
 Formula: C₂₀H₄₆NO₅P
 Molecular Weight: 411.56
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Luviquat mono CP is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4298 mL	12.1489 mL	24.2978 mL
5 mM	0.486 mL	2.4298 mL	4.8596 mL
10 mM	0.243 mL	1.2149 mL	2.4298 mL
50 mM	0.0486 mL	0.243 mL	0.486 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481