

Compound N084-1019

Chemical Properties


CAS No. :

Formula: C₂₄H₄₀O₅

Molecular Weight: 408.579

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N084-1019 is a useful organic compound for research related to life sciences and the catalog number is T131170.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4475 mL	12.2375 mL	24.475 mL
5 mM	0.4895 mL	2.4475 mL	4.895 mL
10 mM	0.2448 mL	1.2238 mL	2.4475 mL
50 mM	0.049 mL	0.2448 mL	0.4895 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481