

## Fumigaclavine C

## Chemical Properties

CAS No. :	62867-47-4
Formula:	C <sub>23</sub> H <sub>30</sub> N <sub>2</sub> O <sub>2</sub>
Molecular Weight:	366.505
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Fumigaclavine C is a useful organic compound for research related to life sciences. The catalog number is T124517 and the CAS number is 62867-47-4.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7284 mL	13.6422 mL	27.2844 mL
5 mM	0.5457 mL	2.7284 mL	5.4569 mL
10 mM	0.2728 mL	1.3642 mL	2.7284 mL
50 mM	0.0546 mL	0.2728 mL	0.5457 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481