

## Dodovisone C

## Chemical Properties


CAS No. :

Formula: C54H58O16

Molecular Weight: 963.042

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Dodovisone C is a useful organic compound for research related to life sciences and the catalog number is T125678.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0384 mL	5.1919 mL	10.3838 mL
5 mM	0.2077 mL	1.0384 mL	2.0768 mL
10 mM	0.1038 mL	0.5192 mL	1.0384 mL
50 mM	0.0208 mL	0.1038 mL	0.2077 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481