

## Dodovisone C

## Chemical Properties

CAS No. :	1616683-52-3
Formula:	C <sub>27</sub> H <sub>30</sub> O <sub>8</sub>
Molecular Weight:	482.52
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Dodovisone C is a useful organic compound for research related to life sciences and the catalog number is T125678.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0725 mL	10.3623 mL	20.7245 mL
5 mM	0.4145 mL	2.0725 mL	4.1449 mL
10 mM	0.2072 mL	1.0362 mL	2.0725 mL
50 mM	0.0414 mL	0.2072 mL	0.4145 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481