

Compound T129334

Chemical Properties

CAS No. :

Formula: C₃₈H₆₈O₁₈

Molecular Weight: 812.94

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	16-(((2R,3R,4R,5S,6R)-6-[(acetyloxy)methyl]-5-(((2S,3R,4S,5R,6R)-3,5-dihydroxy-6-(hydroxymethyl)-4-[(3-hydroxyoctanoyl)oxy]oxan-2-yl)oxy)-3,4-dihydroxyoxan-2-yl)oxy)-2,15-dihydroxyhexadecanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2301 mL	6.1505 mL	12.301 mL
5 mM	0.246 mL	1.2301 mL	2.4602 mL
10 mM	0.123 mL	0.6151 mL	1.2301 mL
50 mM	0.0246 mL	0.123 mL	0.246 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481