


Serratriol

Chemical Properties

CAS No. :	13956-52-0
Formula:	C ₃₀ H ₅₀ O ₃
Molecular Weight:	458.72
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Serratriol is a natural product for research related to life sciences. The catalog number is TN5000 and the CAS number is 13956-52-0.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.180 mL	10.8999 mL	21.7998 mL
5 mM	0.436 mL	2.180 mL	4.360 mL
10 mM	0.218 mL	1.090 mL	2.180 mL
50 mM	0.0436 mL	0.218 mL	0.436 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481