


Bitopertin

Chemical Properties

CAS No. :	845614-11-1
Formula:	C ₂₁ H ₂₀ F ₇ N ₃ O ₄ S
Molecular Weight:	543.46
Storage:	Store at low temperature Powder: -20°C for 3 years In solvent: -80°C for 1 year <i>Actual storage temperature shall be subject to the COA.</i>



Biological Description

Description	Bitopertin (Paliflutine) (RG1678, RO-4917838) is a potent inhibitor of glycine transporter 1 (GlyT1), with K_i of 8.1 nM for human hGlyT1b and IC_{50} of 22-25 nM in Chinese hamster ovary cells.
Targets(IC_{50})	GlyT
In vitro	RG1678 noncompetitively inhibits [³ H]glycine uptake in cells stably expressing hGlyT1b and mGlyT1b, with IC_{50} values of 25 ± 2 nM and 22 ± 5 nM, respectively ($n = 6$) and competitively displaces [³ H]ORG24598 binding with a K_i of 8.1 nM at human hGlyT1b in membranes from Chinese hamster ovary cells. RG1678 has no effect on hGlyT2-mediated [³ H]glycine uptake up to 30 μ M concentration. There is no significant species difference in the pharmacology for RG1678 based on the ability of the compound to displace [³ H]ORG24598. In hippocampal CA1 pyramidal cells, RG1678 enhances NMDA-dependent long-term potentiation (LTP) at 30 nM ($213 \pm 18\%$; $n=7$), 100 nM ($269 \pm 44\%$, $n=13$) but not at 300 nM ($152 \pm 14\%$; $n = 9$)[1].
In vivo	Administration of RG1678 produces a long-lasting (>3h) dose-dependent increase in extracellular glycine levels both in microdialysis experiments conducted in rats and CSF of rats. In mice, RG1678 dose-dependently and significantly attenuates hyperlocomotion induced by the psychostimulant D-amphetamine. RG1678 also prevents the hyper-response to D-amphetamine challenge in rats treated chronically with phencyclidine, an NMDA receptor open-channel blocker[1].
Kinase Assay	Association and dissociation kinetic analysis of [³ H]ORG24598 to hGlyT1 and ratforebrain membranes is performed. [³ H]ORG24598 binding experiments are performed using membranes from CHO cells expressing hGlyT1b and also in membranes from mouse, rat, monkey, and dogforebrains. Saturation isotherms are determined by adding [³ H]ORG24598 to rat, mouse, monkey, and dog forebrain membranes (40 μ g/well) and cell membranes (10 μ g/well) in a total volume of 500 μ L for 3 h at room temperature. Saturation binding experiments are analyzed by an Excel-based curve-fitting program using the Michaelis-Menten equation derived from the equation of a bimolecular reaction and the law of mass action: $B = (B_{max} \times [F]) / (K_d + [F])$, where B is the amount of ligand bound at equilibrium, B_{max} the maximum number of binding sites, [F] the concentration of free ligand, and K_d the ligand dissociation constant. For inhibition experiments, membranes are incubated with 3 nM [³ H]ORG24598 and 10 concentrations of Bitopertin for 1 h at room temperature. Schild

Kinase Assay	increasing concentrations of [3H]ORG24598 (1-300 nM). IC50 values are derived as described above. Ki values are calculated according to the following equation: $K_i = \frac{IC_{50}}{1 + [L]/K_d}$ [1].
Animal Research	Bitopertin (RG1678) is dissolved in H2O with 0.3% Tween 80 (Mice)[1]. Bitopertin (RG1678) is prepared in Polysorbate 80, HEC, Methyl- and Propylparaben pH 6.0 (Rats) [1]. Male NMRI mice (20-30 g) are treated with Bitopertin (0.3, 3, 1, and 10 mg/kg p.o.) or vehicle (p.o.). After 1 min, L-687,414 (50 mg/kg s.c.) or vehicle is given. After 15 min of habituation in the activity chambers, horizontal activity is recorded for 60 min. The time course of Bitopertin effects on L-687,414-induced hyperactivity is also examined; locomotor activity is assessed 2.5, 4.5, and 24 h after administration of Bitopertin (L-687,414 is always given 15 min before the activity procedure). In addition, the effect of subchronic Bitopertin is investigated. Mice receive vehicle or Bitopertin (1 mg/kg p.o.) for 4 consecutive days and L-687,414-induced hyperactivity is evaluated on day 5. Wistar rats receive a 14-day treatment of PCP HC1 (5 mg/kg) or vehicle (NaCl 0.9%, 5 mL/kg i.p.). 24 h following the last injection, rats (6-18 per group) are allowed to individually habituate to the test boxes for 30 min. Rats then received Bitopertin (1, 3, 10 mg/kg p.o.) or vehicle (Polysorbate 80, HEC, Methyl- and Propylparaben pH 6.0; 5 mL/kg p.o.), followed after 1 h by 1 mg/kg D-amphetamine or vehicle i.p. Horizontal activity is recorded directly after the administration of Bitopertin until 120 min after dosing with amphetamine. Data are analyzed by ANOVA supplemented by Fischer's least significant difference post hoc test.

Solubility Information

Solubility	DMSO: 63.75 mg/mL (117.3 mM), Sonication is recommended. H2O: Insoluble, Ethanol: 5 mg/mL (9.2 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8401 mL	9.2003 mL	18.4006 mL
5 mM	0.368 mL	1.8401 mL	3.6801 mL
10 mM	0.184 mL	0.920 mL	1.8401 mL
50 mM	0.0368 mL	0.184 mL	0.368 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Alberati D, et al. Neuropharmacology. 2012, 62(2):1152-1161.
- Bugarski-Kirola D, et al. Eur Neuropsychopharmacol. 2014, 24(7):1024-1036.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481