


(R)-Pyrrolidine-3-carboxylic acid

Chemical Properties

CAS No. :	72580-54-2
Formula:	C ₅ H ₉ NO ₂
Molecular Weight:	115.13
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	(R)-Pyrrolidine-3-carboxylic acid is a proline derivative suitable for biochemical experiments and drug synthesis research.
Targets(IC50)	Amino Acids and Derivatives

Solubility Information

Solubility	H ₂ O: 16.00 mg/mL (138.97 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.6858 mL	43.4292 mL	86.8583 mL
5 mM	1.7372 mL	8.6858 mL	17.3717 mL
10 mM	0.8686 mL	4.3429 mL	8.6858 mL
50 mM	0.1737 mL	0.8686 mL	1.7372 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Luckose F, et al. Effects of amino acid derivatives on physical, mental, and physiological activities. Crit Rev Food Sci Nutr. 2015;55(13):1793-1144.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481