

## Sinococuline

## Chemical Properties

CAS No. :	109351-36-2
Formula:	C <sub>18</sub> H <sub>23</sub> NO <sub>5</sub>
Molecular Weight:	333.38
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sinococuline is a natural product for research related to life sciences. The catalog number is TN6270 and the CAS number is 109351-36-2.
Targets(IC50)	Anti-infection

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9996 mL	14.9979 mL	29.9958 mL
5 mM	0.5999 mL	2.9996 mL	5.9992 mL
10 mM	0.300 mL	1.4998 mL	2.9996 mL
50 mM	0.060 mL	0.300 mL	0.5999 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481