


CBP/p300-IN-19

Chemical Properties

CAS No. : 2592638-13-4
 Formula: C₃₀H₂₇N₃O₃
 Molecular Weight: 477.55
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	CBP/p300-IN-19 is a potent inhibitor of p300/CBP HAT, effective against P300-HAT (IC ₅₀ : 1.4 μM), CBP-HAT (IC ₅₀ : 1.4 μM, 2.2 μM), PCAF (IC ₅₀ : >100 μM), and Myst3 (IC ₅₀ : >100 μM), demonstrating anti-tumor activity.
Targets(IC ₅₀)	Others,Histone Acetyltransferase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.094 mL	10.4701 mL	20.9402 mL
5 mM	0.4188 mL	2.094 mL	4.188 mL
10 mM	0.2094 mL	1.047 mL	2.094 mL
50 mM	0.0419 mL	0.2094 mL	0.4188 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481