


## 1-Indanone

## Chemical Properties

CAS No. :	83-33-0
Formula:	C <sub>9</sub> H <sub>8</sub> O
Molecular Weight:	132.16
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1-Indanone (Indan-1-one) is the intermediate of Indinavir.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 50 mg/mL (378.33 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.5666 mL	37.8329 mL	75.6659 mL
5 mM	1.5133 mL	7.5666 mL	15.1332 mL
10 mM	0.7567 mL	3.7833 mL	7.5666 mL
50 mM	0.1513 mL	0.7567 mL	1.5133 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Fill TP, et al. J Microbiol Biotechnol. 2012 Jun;22(6):832-7.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481