

ANO1-IN-3

Chemical Properties


CAS No. : 2494280-04-3

Formula: C₂₀H₁₇NO₃

Molecular Weight: 319.35

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	ANO1-IN-3 (Compound 3n) is a potent and selective inhibitor of ANO1, with an IC ₅₀ of 1.23 μM, and can induce apoptosis [1].
Targets(IC50)	Apoptosis,Others,Chloride channel

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1314 mL	15.6568 mL	31.3136 mL
5 mM	0.6263 mL	3.1314 mL	6.2627 mL
10 mM	0.3131 mL	1.5657 mL	3.1314 mL
50 mM	0.0626 mL	0.3131 mL	0.6263 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481