


Imitrodast

Chemical Properties

CAS No. :	114686-12-3
Formula:	C ₁₃ H ₁₂ N ₂ O ₂ S
Molecular Weight:	260.31
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Imitrodast is a small molecule thromboxane A ₂ synthase (TXA ₂ synthase) inhibitor for the treatment of immune system disorders, respiratory disorders, and cardiovascular diseases, and can be used in studies of arterial thrombosis, asthma, and coronary vasospasm.
Targets(IC50)	Others,Cytochromes P450,Thrombin

Solubility Information

Solubility	DMSO: 50 mg/mL (192.08 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8416 mL	19.2079 mL	38.4157 mL
5 mM	0.7683 mL	3.8416 mL	7.6831 mL
10 mM	0.3842 mL	1.9208 mL	3.8416 mL
50 mM	0.0768 mL	0.3842 mL	0.7683 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Dogné J M, et al. Thromboxane A₂ inhibition: therapeutic potential in bronchial asthma. American Journal of Respiratory Medicine. 2002;1: 11-17.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481