


## 6-O-Isobutyrylhelenalin, Helenalin isobutyrate

### Chemical Properties

CAS No. :	53658-84-7
Formula:	C <sub>19</sub> H <sub>24</sub> O <sub>5</sub>
Molecular Weight:	332.396
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



### Biological Description

Description	6-O-Isobutyrylhelenalin, Helenalin isobutyrate is a useful organic compound for research related to life sciences. The catalog number is T124696 and the CAS number is 53658-84-7.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0084 mL	15.0421 mL	30.0842 mL
5 mM	0.6017 mL	3.0084 mL	6.0168 mL
10 mM	0.3008 mL	1.5042 mL	3.0084 mL
50 mM	0.0602 mL	0.3008 mL	0.6017 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481