


ONO-8539

Chemical Properties

CAS No. : 459842-29-6
 Formula: C₂₅H₂₈N₂O₅S₂
 Molecular Weight: 500.63
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	ONO-8539 is a prostanooid EP1 receptor antagonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9975 mL	9.9874 mL	19.9748 mL
5 mM	0.3995 mL	1.9975 mL	3.995 mL
10 mM	0.1997 mL	0.9987 mL	1.9975 mL
50 mM	0.0399 mL	0.1997 mL	0.3995 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481