

## Compound T129490

### Chemical Properties


CAS No. :

Formula: C<sub>25</sub>H<sub>26</sub>O<sub>13</sub>

Molecular Weight: 534.47

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description

3-(4-hydroxy-5-methoxy-2-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)phenyl)-2-[(2-oxo-2H-chromen-7-yl)oxy]propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.871 mL	9.3551 mL	18.7101 mL
5 mM	0.3742 mL	1.871 mL	3.742 mL
10 mM	0.1871 mL	0.9355 mL	1.871 mL
50 mM	0.0374 mL	0.1871 mL	0.3742 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481