

BRD-K25923209

## Chemical Properties


CAS No. : 2303501-32-6

Formula: C<sub>29</sub>H<sub>36</sub>N<sub>4</sub>O<sub>3</sub>

Molecular Weight: 488.62

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	BRD-K25923209 is a novel inhibitor of BAF transcriptional repression.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0466 mL	10.2329 mL	20.4658 mL
5 mM	0.4093 mL	2.0466 mL	4.0932 mL
10 mM	0.2047 mL	1.0233 mL	2.0466 mL
50 mM	0.0409 mL	0.2047 mL	0.4093 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481