


## Syringin pentaacetate

## Chemical Properties

CAS No. :	92233-55-1
Formula:	C <sub>27</sub> H <sub>34</sub> O <sub>14</sub>
Molecular Weight:	582.55
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Syringin pentaacetate is a natural product for research related to life sciences. The catalog number is TN5087 and the CAS number is 92233-55-1.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7166 mL	8.583 mL	17.1659 mL
5 mM	0.3433 mL	1.7166 mL	3.4332 mL
10 mM	0.1717 mL	0.8583 mL	1.7166 mL
50 mM	0.0343 mL	0.1717 mL	0.3433 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481