


L 366811

Chemical Properties

CAS No. : 127819-95-8
 Formula: C43H56N8O6
 Molecular Weight: 780.95
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	L 366811 is a potent oxytocin antagonist.
Targets(IC50)	Others,Bradykinin Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2805 mL	6.4025 mL	12.8049 mL
5 mM	0.2561 mL	1.2805 mL	2.561 mL
10 mM	0.128 mL	0.6402 mL	1.2805 mL
50 mM	0.0256 mL	0.128 mL	0.2561 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Pettibone DJ, Clineschmidt BV, Lis EV, Ransom RW, Totaro JA, Young GS, Bock MG, Freidinger RM, Veber DF, Williams PD. Bradykinin agonist activity of a novel, potent oxytocin antagonist. Eur J Pharmacol. 1991 Apr 24;196(3):233-7. PubMed PMID: 1893911.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481