

Lobendazole

Chemical Properties

CAS No. :	6306-71-4
Formula:	C10H11N3O2
Molecular Weight:	205.213
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Lobendazole is a metabolite of thiophanate.
Targets(IC50)	Others,Parasite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8731 mL	24.3653 mL	48.7306 mL
5 mM	0.9746 mL	4.8731 mL	9.7461 mL
10 mM	0.4873 mL	2.4365 mL	4.8731 mL
50 mM	0.0975 mL	0.4873 mL	0.9746 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Delatour P, Besse S. [Ethyl benzimidazole carbamate: teratogenic effect and presence in the milk of cows after administration of thiophanate]. Ann Rech Vet. 1990;21(1):87-92. French. PubMed PMID: 2344141.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481