


## Cholera autoinducer 1

## Chemical Properties

CAS No. : 1004296-82-5  
 Formula: C<sub>13</sub>H<sub>26</sub>O<sub>2</sub>  
 Molecular Weight: 214.34  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Cholera autoinducer 1 is an autoinducer of <i>Vibrio cholerae</i> . Cholera autoinducer 1 influences chemotaxis in <i>C. elegans</i> [1].
Targets(IC50)	Others,Antibacterial

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6655 mL	23.3274 mL	46.6548 mL
5 mM	0.9331 mL	4.6655 mL	9.331 mL
10 mM	0.4665 mL	2.3327 mL	4.6655 mL
50 mM	0.0933 mL	0.4665 mL	0.9331 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481