

Compound T130976

Chemical Properties


CAS No. :

Formula: C34H42O17

Molecular Weight: 722.69

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

(2S,3R,4R,5S,6S)-4,5-bis(acetyloxy)-2-[[[(2S,4S,5S,6R,10S)-2-(hydroxymethyl)-10-[[[(2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]-3,9-dioxatricyclo[4.4.0.0^{2,4}]]dec-7-en-5-yl]oxy]-6-methyloxan-3-yl (2E)-3-phenylprop-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3837 mL	6.9186 mL	13.8372 mL
5 mM	0.2767 mL	1.3837 mL	2.7674 mL
10 mM	0.1384 mL	0.6919 mL	1.3837 mL
50 mM	0.0277 mL	0.1384 mL	0.2767 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481