

Ac-IEVDIDVEH (TFA)

Chemical Properties

CAS No. :

Formula: C50H76F3N11O21

Molecular Weight: 1224.19

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Ac-IEVDIDVEH TFA is a short peptide sequence with Ac at the end.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8169 mL	4.0843 mL	8.1687 mL
5 mM	0.1634 mL	0.8169 mL	1.6337 mL
10 mM	0.0817 mL	0.4084 mL	0.8169 mL
50 mM	0.0163 mL	0.0817 mL	0.1634 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481