


Methanone, 1-piperidinyl(tetrahydro-2H-pyran-4-yl)-

Chemical Properties

CAS No. :	936083-43-1
Formula:	C11H19NO2
Molecular Weight:	197.28
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Methanone, 1-piperidinyl(tetrahydro-2H-pyran-4-yl)-, with CAS No. 936083-43-1, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Methanone, 1-piperidinyl(tetrahydro-2H-pyran-4-yl)- provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0689 mL	25.3447 mL	50.6894 mL
5 mM	1.0138 mL	5.0689 mL	10.1379 mL
10 mM	0.5069 mL	2.5345 mL	5.0689 mL
50 mM	0.1014 mL	0.5069 mL	1.0138 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481