


Annonidine F

Chemical Properties

CAS No. : 820965-40-0
 Formula: C₂₆H₂₈N₂
 Molecular Weight: 368.524
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Annonidine F is a natural product that can be used as a reference standard. The CAS number of Annonidine F is 820965-40-0.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7136 mL	13.5678 mL	27.1356 mL
5 mM	0.5427 mL	2.7136 mL	5.4271 mL
10 mM	0.2714 mL	1.3568 mL	2.7136 mL
50 mM	0.0543 mL	0.2714 mL	0.5427 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481