

Compound T130389

Chemical Properties


CAS No. :

Formula: C₃₂H₄₆O₁₀

Molecular Weight: 590.71

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-{3-[(4E)-6-(acetyloxy)-2-hydroxy-6-methyl-3-oxohept-4-en-2-yl]-6-(carboxymethyl)-2-hydroxy-3a,5a,9b-trimethyl-5-oxo-1H,2H,3H,3aH,4H,5H,5aH,6H,9H,9aH,9bH-cyclopenta[a]naphthalen-7-yl}-2-methylpropanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6929 mL	8.4644 mL	16.9288 mL
5 mM	0.3386 mL	1.6929 mL	3.3858 mL
10 mM	0.1693 mL	0.8464 mL	1.6929 mL
50 mM	0.0339 mL	0.1693 mL	0.3386 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481