

## Compound T128561

## Chemical Properties

CAS No. :

Formula: C<sub>35</sub>H<sub>58</sub>O<sub>9</sub>

Molecular Weight: 622.84

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	3-[3-(2-[[3,4-dihydroxy-5-(hydroxymethyl)oxolan-2-yl]oxy]-5-hydroxy-6-methylhept-6-en-2-yl)-4-hydroxy-6,9a,9b-trimethyl-7-(prop-1-en-2-yl)-dodecahydro-1H-cyclopenta[a]naphthalen-6-yl]propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6055 mL	8.0277 mL	16.0555 mL
5 mM	0.3211 mL	1.6055 mL	3.2111 mL
10 mM	0.1606 mL	0.8028 mL	1.6055 mL
50 mM	0.0321 mL	0.1606 mL	0.3211 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481