

Isobavachin

Chemical Properties


CAS No. : 31524-62-6

Formula: C₂₀H₂₀O₄

Molecular Weight: 324.37

Storage: Keep away from direct sunlight, Store under nitrogen
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Isobavachin can facilitate mouse ES cells differentiating into neuronal cells.
Targets(IC50)	Others, NF-κB, Autophagy, Cytochromes P450, p38 MAPK

Solubility Information

Solubility	DMSO: 59 mg/mL (181.89 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (6.17 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0829 mL	15.4145 mL	30.829 mL
5 mM	0.6166 mL	3.0829 mL	6.1658 mL
10 mM	0.3083 mL	1.5414 mL	3.0829 mL
50 mM	0.0617 mL	0.3083 mL	0.6166 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Wang DY, et al. Promoting effects of isobavachin on neurogenesis of mouse embryonic stem cells were associated with protein prenylation. Acta Pharmacol Sin. 2011 Apr;32(4):425-32.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481