

## Dicresulene diammonium

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Dicresulene diammonium is an impurity found in Policresulen, which is an organic acid possessing hemostatic, antimicrobial, and antiviral properties[1][2].
Targets(IC50)	Others,Antibacterial
In vitro	Dicresulene, a chemical compound utilized in both traditional Chinese and Western medicine, is employed for the treatment of vaginal leukoplakia.

## Reference

Han, Xuping, et al. Traditional Chinese medicine and Western medicine composition for treating leukoplakia of vagina and preparation method thereof. CN108283688A.

Deng-wei Wu, et al. Policresulen, a novel NS2B/NS3 protease inhibitor, effectively inhibits the replication of DENV2 virus in BHK-21 cells. Acta Pharmacol Sin. 2015 Sep;36(9):1126-36.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481