


Tinoridine

Chemical Properties

CAS No. :	24237-54-5
Formula:	C ₁₇ H ₂₀ N ₂ O ₂ S
Molecular Weight:	316.42
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tinoridine is a quinoline derivative analog that acts by inhibiting the production of cyclooxygenase-2 (COX-2). It has anti-inflammatory, anti-allergic and anticoagulant activities.
Targets(IC50)	COX,Cytochromes P450

Solubility Information

Solubility	DMSO: 50 mg/mL (158.02 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1604 mL	15.8018 mL	31.6036 mL
5 mM	0.6321 mL	3.1604 mL	6.3207 mL
10 mM	0.316 mL	1.5802 mL	3.1604 mL
50 mM	0.0632 mL	0.316 mL	0.6321 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481