

Delcosine

Chemical Properties


CAS No. :

Formula: C₂₄H₃₉N₇O

Molecular Weight: 453.576

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Delcosine is a natural product that can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2047 mL	11.0234 mL	22.0468 mL
5 mM	0.4409 mL	2.2047 mL	4.4094 mL
10 mM	0.2205 mL	1.1023 mL	2.2047 mL
50 mM	0.0441 mL	0.2205 mL	0.4409 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481