


XE 820

Chemical Properties

CAS No. : 71791-90-7
 Formula: C₃₅H₄₀N₄O₆
 Molecular Weight: 612.72
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	XE 820 is a biochemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6321 mL	8.1603 mL	16.3207 mL
5 mM	0.3264 mL	1.6321 mL	3.2641 mL
10 mM	0.1632 mL	0.816 mL	1.6321 mL
50 mM	0.0326 mL	0.1632 mL	0.3264 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481