


IDN-1965

Chemical Properties

CAS No. : 204919-14-2
 Formula: C₂₁H₂₆N₃O₅
 Molecular Weight: 419.45
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	IDN-1965 is a potent, broad-spectrum, irreversible caspase inhibitor. IDN-1965 is markedly effective at inhibiting Fas-mediated apoptosis by multiple routes of administration. The therapeutic potential of caspase inhibitors appears promising for the treatment of apoptosis-mediated liver injury based on potency and postinsult efficacy.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3841 mL	11.9204 mL	23.8407 mL
5 mM	0.4768 mL	2.3841 mL	4.7681 mL
10 mM	0.2384 mL	1.192 mL	2.3841 mL
50 mM	0.0477 mL	0.2384 mL	0.4768 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481