


## YycG inhibitor 16

## Chemical Properties

CAS No. : 506429-34-1  
 Formula: C<sub>30</sub>H<sub>33</sub>N<sub>3</sub>O<sub>5</sub>  
 Molecular Weight: 547.67  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	YycG inhibitor 16 is a novel YycG inhibitor, inhibiting planktonic cells of <i>E. faecalis</i> , including the vancomycin- or linezolid-resistant strains.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8259 mL	9.1296 mL	18.2592 mL
5 mM	0.3652 mL	1.8259 mL	3.6518 mL
10 mM	0.1826 mL	0.913 mL	1.8259 mL
50 mM	0.0365 mL	0.1826 mL	0.3652 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481