


Pz-1

Chemical Properties

CAS No. : 1800505-64-9
 Formula: C₂₆H₂₆N₆O₂
 Molecular Weight: 454.52
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pz-1 is an inhibitor of VEGFR2 and RET (rearranged during transfection) tyrosine kinase that blocks the blood supply required for RET-stimulated growth.
Targets(IC50)	c-RET, VEGFR
In vitro	Pz-1 (1.0nM) strongly inhibited phosphorylation of all tested RET oncoproteins in cell-based assays with IC50s of less than 1 nM for both VEGFR2 and RET wild type kinases[1].
In vivo	Pz-1 (1.0mg/kg) abrogated the formation of tumors induced by RET-mutant fibroblasts and blocked the phosphorylation of both RET and VEGFR2 in tumor tissue. Pz-1 featured no detectable toxicity at concentrations of up to 100.0mg/kg[1].

Solubility Information

Solubility	DMSO: 45 mg/mL (99.01 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (4.4 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2001 mL	11.0006 mL	22.0012 mL
5 mM	0.440 mL	2.2001 mL	4.4002 mL
10 mM	0.220 mL	1.1001 mL	2.2001 mL
50 mM	0.044 mL	0.220 mL	0.440 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Frett B, et al. Fragment-Based Discovery of a Dual pan-RET/VEGFR2 Kinase Inhibitor Optimized for Single-Agent Polypharmacology. *Angew Chem Int Ed Engl.* 2015 Jul 20;54(30):8717-21.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481