

3-(3,4-dihydroxyphenyl)-2-oxopropyl (2E)-3-(3-hydroxy-4-methoxyphenyl)prop-2-enoate

Chemical Properties


CAS No. :

Formula: C₁₉H₁₈O₇

Molecular Weight: 358.35

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

3-(3,4-dihydroxyphenyl)-2-oxopropyl (2E)-3-(3-hydroxy-4-methoxyphenyl)prop-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7906 mL	13.9528 mL	27.9057 mL
5 mM	0.5581 mL	2.7906 mL	5.5811 mL
10 mM	0.2791 mL	1.3953 mL	2.7906 mL
50 mM	0.0558 mL	0.2791 mL	0.5581 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481