

SAPO 50

Chemical Properties


CAS No. : 193804-76-1

Formula: C75H118O39

Molecular Weight: 1643.73

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	SAPO 50 is a natural product that can be used as a reference standard. The CAS number of SAPO 50 is 193804-76-1.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6084 mL	3.0419 mL	6.0837 mL
5 mM	0.1217 mL	0.6084 mL	1.2167 mL
10 mM	0.0608 mL	0.3042 mL	0.6084 mL
50 mM	0.0122 mL	0.0608 mL	0.1217 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481