


LY063518

Chemical Properties

CAS No. : 182633-64-3
 Formula: C₁₂H₁₁ClN₄S
 Molecular Weight: 278.76
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	LY063518 is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5873 mL	17.9366 mL	35.8732 mL
5 mM	0.7175 mL	3.5873 mL	7.1746 mL
10 mM	0.3587 mL	1.7937 mL	3.5873 mL
50 mM	0.0717 mL	0.3587 mL	0.7175 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481