

## 3'-Azido-3'-deoxycytidine

## Chemical Properties


CAS No. : 70580-87-9

Formula: C<sub>9</sub>H<sub>12</sub>N<sub>6</sub>O<sub>4</sub>

Molecular Weight: 268.23

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Nucleoside Derivatives -Azido-nucleosides, 3'-Modified nucleosides; Drugs and Inhibitors; Ribonucleotide reductase inhibitor
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7281 mL	18.6407 mL	37.2814 mL
5 mM	0.7456 mL	3.7281 mL	7.4563 mL
10 mM	0.3728 mL	1.8641 mL	3.7281 mL
50 mM	0.0746 mL	0.3728 mL	0.7456 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481