

24-Deacetylalisol O

Chemical Properties


CAS No. :

Formula: C₃₀H₄₆O₄

Molecular Weight: 470.694

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	24-Deacetylalisol O is a useful organic compound for research related to life sciences and the catalog number is T126460.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1245 mL	10.6227 mL	21.2454 mL
5 mM	0.4249 mL	2.1245 mL	4.2491 mL
10 mM	0.2125 mL	1.0623 mL	2.1245 mL
50 mM	0.0425 mL	0.2125 mL	0.4249 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481