


2-Hydroxy-6,7-dichloroquinoxaline

Chemical Properties

CAS No. :	78470-95-8
Formula:	C ₈ H ₄ Cl ₂ N ₂ O
Molecular Weight:	215.036
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Hydroxy-6,7-dichloroquinoxaline, with CAS No. 78470-95-8, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2-Hydroxy-6,7-dichloroquinoxaline provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6503 mL	23.2515 mL	46.503 mL
5 mM	0.9301 mL	4.6503 mL	9.3006 mL
10 mM	0.465 mL	2.3251 mL	4.6503 mL
50 mM	0.093 mL	0.465 mL	0.9301 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481