

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0343 mL	10.1717 mL	20.3434 mL
5 mM	0.4069 mL	2.0343 mL	4.0687 mL
10 mM	0.2034 mL	1.0172 mL	2.0343 mL
50 mM	0.0407 mL	0.2034 mL	0.4069 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Pehrson R, et al. AZD7624, an Inhaled p38 Inhibitor, Demonstrates Local Lung Inhibition of LPS-Induced TNF α with Minimal Systemic Exposure. J Pharmacol Exp Ther. 2018 Jun;365(3):567-572.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481