

Fmoc-Dab(Boc)-OH

Chemical Properties

CAS No. :	125238-99-5
Formula:	C ₂₄ H ₂₈ N ₂ O ₆
Molecular Weight:	440.489
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Fmoc-Dab(Boc)-OH is a chemical intermediate used for chemical synthesis.
Targets(IC50)	Others,Amino Acids and Derivatives

Solubility Information

Solubility	DMSO: Soluble, H ₂ O: Insoluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2702 mL	11.351 mL	22.702 mL
5 mM	0.454 mL	2.2702 mL	4.5404 mL
10 mM	0.227 mL	1.1351 mL	2.2702 mL
50 mM	0.0454 mL	0.227 mL	0.454 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Yamada K, Urakawa H, Oku H, Katakai R. Facile synthesis of Nalpha-protected-L-alpha,gamma-diaminobutyric acids mediated by polymer-supported hypervalent iodine reagent in water. J Pept Res. 2004 Aug;64(2):43-50. PubMed PMID: 15251030.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481