


## 3-AQC

## Chemical Properties

CAS No. :	201216-42-4
Formula:	C <sub>20</sub> H <sub>21</sub> N <sub>5</sub> O <sub>4</sub>
Molecular Weight:	395.41
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-AQC is a competitive antagonist of 5-HT <sub>3</sub> .
Targets(IC <sub>50</sub> )	5-HT Receptor

## Solubility Information

Solubility	DMSO: 11 mg/mL (27.82 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.529 mL	12.6451 mL	25.2902 mL
5 mM	0.5058 mL	2.529 mL	5.058 mL
10 mM	0.2529 mL	1.2645 mL	2.529 mL
50 mM	0.0506 mL	0.2529 mL	0.5058 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481