

Isoretronecanol

Chemical Properties


CAS No. :

Formula: C₈H₁₅NO

Molecular Weight: 141.214

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Isoretronecanol is a useful organic compound for research related to life sciences and the catalog number is T124711.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.0817 mL	35.4083 mL	70.8165 mL
5 mM	1.4163 mL	7.0817 mL	14.1633 mL
10 mM	0.7082 mL	3.5408 mL	7.0817 mL
50 mM	0.1416 mL	0.7082 mL	1.4163 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481